Continue



Temparin max instructions

English Please consult a doctor prior to use in pregnant/nursing women, and children. Caution: To avoid infection and/or serious injury, this product should not be used if there is throbbing pain, ulceration or swelling in the affected area. If any of these symptoms exist, consult your dentist immediately. This product should not be placed on open wounds. Do not swallow. Do not use if allergic to Eugenol or Sesame Oil. READ ALL OF THE DIRECTIONS CAREFULLY BEFORE USE. Wash and dry your hands. Use only what is needed to repair the dental emergency, and save the rest for a future application. Keep out of reach of children. Do not get into eyes Tooth Saver: If tooth is dislodged, collect the tooth and do not try to clean it, there could be important soft tissue that will help in the healing process. Fill Tooth Saver container with milk or your own saliva. Place broken or whole tooth in container making certain that milk or saliva is completely covering tooth. Contact your dentist immediately. Chances of successful replantation are best if you see the dentist within 1 hour of incident. Temparin Max : This is a temporary remedy meant to replace a lost filling, cap, crown or inlay until you have a chance to visit your dentist and should only be removed by a dental professional. Visit dentist within 48 hours of use. If irritation persists, inflammation develops, or if fever and infection develop stop use and ask a dentist or doctor. Directions to replace a lost filling: Wash the cavity area with warm water to remove food particles. Do not dry the area. Saturate the tip of a fresh cotton swab applicator with liquid from the Eugenol vial and apply to tooth. The Eugenol will help the filling material set up and adhere to the tooth enamel. Avoid touching surrounding gum tissue, cheek, or tongue with wet applicator tip. Do not leave the applicator in mouth. Cut with spade tool or pinch a small amount of Temparin Max from the vial, and roll into a ball between fingers. Slightly overfill the cavity with Temparin Max and tamp into place with applicator tool tip. Close your teeth together, as if you were biting down on a piece of paper, and gently chew or grind on the replaced filling a couple times to create a comfortable bite. Use a fresh moist cotton swab to remove excess Temparin Max from around the area. Temparin Max will set softly in your mouth within a few minutes, but avoid chewing on the repair for 2-3 hours to allow the material to fully set. Seal any remaining Temparin Max back into the vial to avoid it hardening over time. Directions to secure loose crowns, caps, or inlays: Try the crown, cap, or inlay first to see if you can fit it accurately back onto the tooth surface. If you are not able to do so, or if remaining tooth is cracked, then do not proceed further. If the crown has a post attached to it or if there is tooth surface remaining inside the crown, Temparin Max may not adhere. Carefully rinse the tooth area and the inside of the crown, cap, or inlay with warm water, and do not dry. Cut with spade tool or pinch a small amount of Temparin Max from the vial, roll into a ball between fingers, and place evenly inside the crown, cap, or inlay. Place the crown, cap, or inlay on the tooth and bite down, applying pressure to push it into position. Temparin Max will begin to harden within a few minutes, but avoid chewing on the repair for at least 2-3 hours to allow the material to set. Use a toothbrush to gently clear any excess material from around the crown, cap, or inlay. Seal any remaining Temparin Max back into the vial to avoid it hardening over time. Please consult a doctor prior to use in pregnant/nursing women, and children. Ingredients for Temparin Max: Zinc Oxide Powder, Calcium Sulphate, Petroleum Jelly, Potassium alum, Paraffin, Aluminum Sulphate, Aluminum Phosphate, Menthol Crystals, Eugenol, Yellow Iron Oxide Ingredients for Eugenol: Eugenol, Sesame Oil Caution: To avoid infection and/or serious injury, this product should not be used if there is throbbing pain, ulceration or swelling in the affected area. If any of these symptoms exist, consult your dentist immediately. This product should not be placed on open wounds. Do not swallow. Do not use if allergic to Eugenol or Sesame Oil. The current Corona-virus crisis has a lot of us doing things we never thought we would and DIY dentistry is one of them! I had a filling drop out (whilst eating a rather succulent lamb cutlet we had cooked on the BBQ) and whilst I wasn't in any pain it was going to be some time before I could see a dentist. After reading a few reviews we thought we'd try this kit on the basis it was unlikely to do any harm. If there was a disappointment it was in the length of time it took to arrive (more than10 days) but once here Mrs M gamely said she'd give it a go. It was quite easy despite her trepidation and the temporary filling is still in situ a week later. Let's hope it last a good deal longer ... Sceptical as I've used stuff like this before and they didn't work. However this has the bonding agent which others didn't. I had a filling fall out and needed something as the dentist isn't open. The bonding agent smells (and tastes) like a really strong clove oil. It makes you really salivate so prepare for drooling! I don't think this would work for a chipped tooth as I don't think there would be enough grip to stay put. But as I've got the hole drilled by a dentist so there is a good area to stuff the putty into if that makes sense. Just put the bonding agent in the hole and stuff the putty in with the tool. Bite and grind to smooth. I did mine last thing at night after a good brush and clean out of the hole and next morning it was set. It's been in for a couple of weeks now. Able to eat on it. Brush and it's not moved yet. I understand it's a temporary fix, but before I could not bite on that side due to pressure pain (no pain without biting on it) and now I can with no pain. For a temp fix in these times it's a great product and does the job Ideal emergency kit at this time of health crisis The kit was easy to use and has stayed in place for a couple of weeks now with no issues. However the eugenol you dab the tooth with before applying the putty is very pungent and leaves a strong smell and taste in the mouth which is unpleasant whilst you're applying the filling. If you have a drink of water straight after it clears, so no real issue. You need much less putty than you'd think, so there are plenty of uses in this kit. Great value and has saved all tooth ache and sensitivity! This literally is a life saver when you have a broken tooth or fallen out filling and are in a lot of pain. It's a temporary solution but lasts long enough until you see your dentist. Never knew that it was existed! I had to buy it during lockdown where I had a broken filling. Its doing a really good job! I was told that it lasts for a couple of days, but on my tooth last for a week! Then I gad to apply it again! Bought 2 and i used 1 1/2 of them during a 2months period! Totally recommend it!!! Now that I know it exists such an amazing thing, I will make sure that I always have one in case of emergency I had read lots of reviews about this item, and tried it this morning. So far so good! When it says a small amount, it really does mean tiny. Impressed with how easy it was. Will definitely keep as part of the household emergency supply! With no dentists working under Covid 19 lockdown it was a great help in fixing the sharp remains of a broken tooth. We use cookies on our website to give you the most relevant experience by remembering your preferences. By clicking "Accept", you consent to the use of ALL the cookies. View our policyManage consent About the author: Dr David Chen, DDS Hello, I'm Dr Chen and I'm an actively practicing dentist in Long Island City, NY. I graduated from Columbia University College of Dental Medicine in 2016 but prior to going to dental school I was already working in the dental field. It's been more than a decade since I first got to know dentistry and let me tell you, time flies by quickly. Since then I've developed a fondness for writing, which is how this all got started! Association Memberships: American Dental Association New York State Dental Association Queens Count Dental Society Medical Disclaimer: This blog is purely meant for information purposes and should not be used as medical advice. Each situation in your mouth is unique and complex. It is not possible to give advice nor diagnose any oral conditions based on text nor virtual consultations. The best thing to do is to go in person to see your dentist for an examination and consultation so that you can receive the best care possible. The purpose of all of this oral health information is to encourage you to see your dentist and to inform you of what you may expect during your visit. Due to the unfortunate nature of dentistry, there isn't really any true home remedies that will get rid of dental problems. Roughly 99.99% of them require in-person intervention by a healthcare professional. Hint: That is the reason why you can't eliminate seeing dentists in your life! English It is recommended to visit your dentist within 48 hours of application. Please consult a doctor prior to use in pregnant/nursing women, and children. Caution: To avoid infection and/or serious injury, this product should not be used if there is throbbing pain, ulceration or swelling in the affected area. If any of these symptoms exist, consult your dentist immediately. This product should not be placed on open wounds. Do not use if allergic to Eugenol. READ ALL OF THE DIRECTIONS CAREFULLY BEFORE USE! Wash and dry your hands. Do not get unused Temparin Max wet. Use only what is needed to repair the dental emergency, and save the rest for a future application. KEEP OUT OF REACH OF CHILDREN. Directions to replace a lost filling: Wash the cavity area with warm water and do not dry. Cut with spade tool or pinch a small amount of Temparin Max from the vial, and roll into a ball between fingers. Slightly overfill the cavity with Temparin Max and tap into place with applicator tool tip. Wet your finger and dab down on the Temparin Max while it is still soft to create a smooth surface. Close your teeth together, as if you were biting down on a piece of paper, and gently chew or grind on the replace filing a couple times to create a comfortable bite. Use a moist cotton swab to remove excess Temparin Max from around the area. Temparin Max will set softly in your mouth within a few minutes, but avoid chewing on the repair for 2-3 hours to allow the material to fully set. Seal any remaining Temparin Max back inside the vial to avoid it hardening over time. Directions to secure loose crowns, caps or inlays: Try the crown, cap or inlay first to see if you can fit it accurately back onto the tooth surface. If you are not able to do so, or if the remaining tooth is cracked, then do not proceed further. If the crown has a post attached to it or if there is tooth surface remaining inside the crown, Temparin Max may not adhere. Carefully rinse the tooth area and the inside of the crown, cap or inlay with warm water, and do not dry. Cut with spade tool or pinch a small amount of Temparin Max from the vial, roll into a ball between fingers and place evenly inside the crown, cap or inlay. Place the crown, cap or inlay on the tooth and bite down, applying pressure to push into position. Temparin Max will begin to harden within a few minutes, but avoid chewing on the repair for at least 2-3 hours to allow the material to set. Use a toothbrush to gently clear any excess material from around the crown, cap or inlay. Seal any remaining Temparin Max back inside the vial to avoid it hardening over time. Please consult a doctor prior to use in pregnant/nursing women, and children. Caution: To avoid infection and/or serious injury, this product should not be used if there is throbbing pain, ulceration or swelling in the affected area. If any of these symptoms exist, consult your dentist immediately. This product should not be placed on open wounds. Do not use if allergic to Eugenol. READ ALL OF THE DIRECTIONS CAREFULLY BEFORE USE! Wash and dry your hands. Do not get unused Temparin Max wet. Use only what is needed to repair the dental emergency, and save the rest for a future application. KEEP OUT OF REACH OF CHILDREN. amazing product simple easy to use if you follow instructions! so after reading reviews i thought thought this wouldnt work but it was cheap so thought id try it out ... i had a hole in one tooth and lost filling from another which was causing me the worst tooth ache ive ever had !!! .. what you need to remember is this isnt the same thing they use to fill a tooth at the dentist so it takes a while to set ! .. i brushed my teeth before then rolled up a small bit of the filling i used my finger to apply rather than the plastic tool as i wanted to press it in and make sure it was covered then just added another ball and pressed that it the i ran my finger under warm water and rubbed it over it to help set it in place .. i found when i closed my mouth it was a good idea to press down my top jaw to bottm to make sure there wasnt too much in and i would be able to chew ok .. after doing this i opened my mouth and with a toothbrush scrapped the little bits that had come loose of my tongue and sides of my mouth .. i applied this at night before i went to bed DO NOT TRY AND DO IT IN THE DAY TIME IT WONT WORK !! .. it needs a long time to set in place .. i went to bed at 11 and was up at 6am and it had set perfectly . the next day i tried to eat soft foods but the day after that i was eating pizza and crisps and ive had no problems at all .. the toothache i had was gone by the next morning as well i couldnt belie how well it worked ! .. the tooth with the hole ive had to re do a number of times as the hole isnt round its on an edge so it doesnt pack in as well and does crumble out if i crunch something very hard but thats to be expected ! but the tooth i filled with a hole in the middle of the tooth has never fallen out set hard and has been in place now for over a month ! .. to me an amazing product if you take the time to do it right and let it set . My dentist was adamant one of my teeth could not be filled and started telling me how I need a private appointment and thousands of pounds to fix a few of my teeth and could sense a real conflict of interest. I bought this product, shoved some in my back tooth and two weeks later its still there doing just fine. Life Saver!d This stuff is just amazing, just wash mouth out, dry the area that needs filling, get some of the filler out roll into a ball and fill the hole! It gets hard really fast and can last weeks if needed! I always have one of these in the house at all times, as you never know when you will need it This filling lasted 10 days max- went to dentist-their filling lasted 3 weeks- so I'm putting this filling Den Tek back in till I can get the filling done again, although not till after Christmas now,it's good stuff, The method is to roll a small amount on to the affected tooth which must be cleaned &dried , then nothing to eat of that side for at least 2 hours-liquids too or it will come away As the others reviewers have said this is only a temporary fix until you can get to the dentist but great in an emergency. It does stay relatively soft and although this means you have to top it up from time to time at least your dentist won't have to drill it out! On but if advice, it's tempting to put too much in the gap but it's better to use small amounts. We bought this as a temporary measure until the dentist could fit us in. Used for one small hole in a molar. Sometimes hard to make it stay initially but has lasted really well and seems to be preventing the pain. It's a tiny pot (although it's plenty) DenTek Temparin Max: Lots of reviews say this stuff is soft when set, but that's a bit misleading. I'd say it's about the same as plaster of paris – nowhere near as hard as your tooth enamel, but not exactly soft either. After all, it does need to be easy enough for your dentist to remove when (if) you go for proper treatment. Really easy to use. Just cut a bit off and work it between thumb and forefinger to soften it up, roll it into a ball, then squeeze it into the cavity with your finger tip, where it begins to harden on contact with water. It doesn't stick to your tooth as such, just moulds itself to irregularities in the cavity to stay in place. I apply it last thing at night and when I wake up in the morning it's set hard and stays in place for a week or two. You only get a small amount – a stick about 15mm long by 7mm in diameter, but that's enough to patch up the very big hole in my back tooth at least 6 times I would say. The only truly bad thing I can say about this stuff is it's giving me an excuse to put off going to the dentist! DenTek First Aid Kit – Emergency Tooth Repair View product We use cookies on our website to give you the most relevant experience by remembering your preferences. By clicking "Accept", you consent to the use of ALL the cookies. View our policyManage consent Privacy Policy If you no longer want to be a part of Club Publix, you can opt out of Club Publix here. At Publix Super Markets, Inc. (collectively with or its affiliates, brands and subsidiaries, "Publix", "we", "our", "us" or any derivatives thereof), we are committed to protecting your privacy. We developed this Privacy Policy (the "Policy") to inform you about how we treat the information that we collect about you in association with our websites, mobile applications, and other online and offline services (collectively, the "Services"). This Policy also describes your choices about the collection and use of your information. Please read this policy carefully before you start to use our Services. By using the Services, you consent to our posted Terms of Use (the "Terms") together with this Policy, which the Terms incorporate by reference. If you do not agree to our Terms with this Policy, or if you violate them in any way, your right to access or use the Services is terminated. PLEASE SEE OUR TERMS REGARDING YOUR RIGHTS IN ANY DISPUTE INVOLVING OUR SERVICES, INCLUDING A MANDATORY ARBITRATION CLAUSE, AND A WAIVER TO A JURY TRIAL OR CLASS ACTION, THAT LIMITS YOUR RELIEF. NOTE TO PHARMACY CUSTOMERS: Notwithstanding the disclosures in this Policy, all pharmacy activities that involve Protected Health Information ("PHI") as defined under the Health Insurance Portability & Accountability Act and related laws and regulations (collectively referred to as "HIPAA") are subject to our posted HIPAA Notice of Privacy Practices. Publix reserves the right, at any time, to modify this Policy. By continuing to use the Services following such modifications, you agree to be bound by such modifications. We will post our most current Policy on this page with the "last updated" date, and you should periodically visit this page to review it. 1. How We Collect Information We may collect information about you by various means, including: Directly from you, both online (e.g., our website or app) and offline (e.g., in store or by phone); From other third-party sources and social media platforms that you may use to engage with us, as well as from partners or third parties we work with; and By combining information from different sources, including online and offline data. 2. Categories of Information That We Collect a. General Categories While operating our business, we collect a range of personal information. That information includes: Identifiers, such as your name, physical address, email address, social media handle, Club Publix and GreenWise Market identifiers, phone number, IP address, device identifiers, signature, government identifiers collected for certain sales and ID verification, and identifiers associated with cookies, pixels or other similar technologies; Financial information, such as payment card information; Health and insurance information, such as prescription information if you are a customer of our pharmacies; Demographics, such as your age and gender; Transactional information, such as records of your purchases and orders, products you like, coupons or offers you select, and content you access, download, or view via the Services; Internet or other electronic network activity, such as information about how you use, access, or interact with the Services or our emails, including but not limited to information about how you use your mouse, keyboard, and other input devices when accessing the Services; information about your activities on other online services; and information about your browser and device; Location information, such as city and state you provide, or location information derived from your device or browser, subject to the settings on your device or browser; Records of communications with our customer care teams or digital assistants for training, recordkeeping, fulfillment, and quality assurance purposes; Video recordings in public areas for security, safety, asset protection and other operational purposes, such as measuring traffic patterns and protecting and improving our business; Survey responses; and Inferences regarding your preferences or other characteristics based on the information we collect or obtain from other sources. Please read below for more information about the information we collect. b. Information That We Collect from You You may be able to use some of our Services or conduct in-store business with us without providing any personal information to us. However, in order to take advantage of certain features, offerings, or other functions available in stores or on our Services, you may be asked or required to provide personal information. For example, depending on the feature or service you wish to access, we may ask you to register an account with us, provide identification, or otherwise collect information from you such as your name, address, e-mail address, telephone number, credit card information, household information, demographic information, and/or information about your interests. You can always choose to not provide us with the requested information; however, if you choose to not provide such information, you may not be able to engage in certain in store transactions or access or use certain features or functions of the Services. When you shop in our stores or on our website, we may collect information about the products or services you buy when you use a payment card, the Publix App or enter a number at the register. To the extent that you provide us with a mobile telephone number you hereby: represent that you are the accountholder and subscriber for the mobile telephone number that you provide, or that you have the accountholder and subscriber's permission to receive text messages and/or phone calls at that telephone number; agree to receive, and authorize us to deliver or cause to be delivered to you, calls and SMS/text messages for any purpose, including advertising, using an automatic telephone dialing system or an artificial or pre-recorded voice; agree to notify us promptly if the mobile telephone number you have provided changes; agree to indemnify us in full for all claims, expenses, costs, and damages related to, arising under, or caused by, in whole or in part, your failure to notify us promptly if the phone number you have provided changes, including but not limited to, all claims, expenses, costs, and damages related to or arising under the Telephone Consumer Protection Act; and understand and agree that your agreement to these terms is not required as a condition of purchasing any property, goods, or services and that, based on your plan with your mobile carrier, normal usage rates and fees, including data fees and text message rates, may apply to text messages and/or phone calls delivered by Publix or on Publix's behalf. c. Information We Collect from Other Sources We may also receive information about you from other sources, including third parties that help us: update, expand, and analyze our records; identify new customers; collect information on our behalf; and prevent or detect fraud. We may also receive information about you from social media platforms websites and third-party applications, including but not limited to when you interact with us on those platforms, access our social media content, or use your social media login credentials to access our Services. The information we receive from these sources may include, but is not limited to, identifiers, demographics, education information, address information, Internet or other electronic network activity, and inferences about your preferences or other characteristics. The information we will be able to receive will be determined by the privacy settings, policies, and/or procedures of the applicable third-party platform, website or service. d. Information We Collect Automatically When you use the Services, we or third-parties may use cookies and similar technologies to collect information about how you access, use, or otherwise interact with the Services, our emails, and online advertising. This information is collected for various reasons, including to help ensure that the Services continue to work correctly, to improve the Services, to support marketing campaigns (including targeted advertising), to understand our customers and how they use the Services, and to offer you a more personalized experience. The information collected may include technical information such as your IP address, the type of device or browser that you use to access the Services, the operating system on your device, the website or other online service that directed you to our Services, the URLs of the websites you visited before visiting our site, numeric identifiers for your mobile device, identifiers for advertising, identifiers generated from secure hash algorithms, session times and durations, how you use your mouse, keyboard, and other input devices to interact with the Services, and the content that you view on the Services. For more information about how this information is collected, please see our Cookies and Other Tracking Technologies section below. e. Information Collected from Mobile Devices When you access our Services via a mobile device, Publix or its third-party service providers may collect various information, including the type of mobile device you use, the temporary or persistent device identifier associated with your device, the IP address of your mobile device, the device's operating system, the browser you use to access the Services, information about how you use the Services, and your device's location as described below. You may be able to use our mobile applications to scan barcodes or take images of products or prescription bottles in order to create shopping lists or conduct transactions. To do so, our mobile applications may request permission to access your camera and photo files. Please note that some features of our mobile applications and/or Services may not be usable or available to you if you do not grant us such permission. f. Location Information Our mobile applications and websites may also collect precise information about the location of your device after you permit us to obtain such information from technologies such as Bluetooth beacons, GPS, Wi-Fi signals, or cell tower signals. We may also use technologies to determine if you are shopping in one of our stores by relying on signals transmitted by your mobile devices. We collect this information to provide you with services, content and offers based on your location. For example, we may use our mobile applications to provide you with the location of nearby stores and deals available in those stores. You may be able to enable our Services to notify us when you enter one of our stores so that we can prepare any orders you may have or to let you know about the events and deals available in the store. To opt out of the collection of precise location or presence information, you can adjust the permissions on your mobile device or uninstall our mobile applications. Publix has no control over your device's settings. And adjusting these settings does not completely eliminate the ability of others to develop location information about your device. g. Protected Health Information We understand that the privacy and security of your health information is important to you. As noted above, when you shop or interact with our pharmacies, we have procedures in place to protect your health information in accordance with our Publix Pharmacy Notice of Privacy Practices. h. Information You Give Us About Others You may give us information about other people, such as the name and address of a gift recipient, or the name and contact information of someone who will pick up items for you at a store. We will only use this information for the specific reason you provide it. 3. Why Do We Collect Your Information We use the information we collect for purposes allowed by applicable laws, including to: fulfill your requests for our products, programs, or services; complete and manage purchases, orders, refunds, returns and exchanges; respond to your inquiries; offer you other products, programs, or services that we believe may be of interest to you; achieve the purposes for which you provided the information; conduct surveys, research, analysis or market studies; operate, evaluate and enhance our products, services, or business operations; prevent or detect fraud; deliver promotional and marketing communications, special offers, or advertisements you via physical mail and electronically that may be of interest to you, including the delivery of targeted advertisements that are based on your previous online activity on the Services and other third-party services (for example, if you view a page on our website about a particular product, you may receive an online advertisement for that product or a related product on a different area of our Services or on other websites on the Internet you visit); create and deliver personalized content, offers, services, and advertising, including personalized content based on the precise location information we collect from mobile devices or the location that we infer from your IP address or other information; maintain and administer your Club Publix or GreenWise Markets accounts or other accounts or memberships you may have with us; conduct promotional activities, including contests; communicate with you regarding your accounts, purchases, or other issues; protect the rights or safety of you, us or others; and comply with and enforce applicable legal requirements, contractual obligations (including our contractual obligations with you, such as through our Club Publix and GreenWise Markets programs, if applicable, or digital properties' terms of use), relevant industry standards, and our policies. We may combine the information we collect from you with information that we collect from other sources. We may use and disclose non-personal, non-individual statistics or demographic information in aggregate form without restriction. 4. How We May Disclose Your Information We may share any of the categories of personal information we collect in the following circumstances: We may share information with business partners and advertisers for our marketing and advertising purposes, but we do not monetize your personal information by selling or renting it for monetary compensation to outside parties; We may share your personal information to our service providers, who use that information to perform technology, business or other professional functions on our behalf, including but not limited to, hosting our Services, delivering transactional communications, fulfilling orders, and processing payments. We only provide such vendors with information so they can perform their required functions on our behalf; We may share anonymized information with third parties once it is altered, combined, or otherwise processed so that it can't reasonably be identified with an individual person; We may disclose the information to third parties when we are required to do so by any law, regulation, court order, subpoena, law enforcement, or other regulatory requirement; We may disclose the information if we believe that such disclosure is necessary or appropriate to comply with applicable law; operate the Services; protect the rights, property or safety of Publix, other users of the Services, partners, affiliates, sponsors, providers, licensors, merchants, or others; prevent harm or loss; or in connection with an investigation of suspected or actual unlawful activity; We may share information with our related subsidiaries and affiliates, subject to the terms of this Policy; We may disclose the information in association with any potential or actual sale or other transfer of some or all of our business assets; and We may disclose information subject to your consent or as necessary to fulfill the purposes for which the information was collected. Where appropriate, we will limit sharing of your information in accordance with the choices you have provided us and applicable law. 5. Public Information If you post comments, pictures, or other content via public forums that may be available on the Services, that information may be publicly available. The information may be searchable by search engines and copied or republished by third parties. We have no control over those third parties. Please think carefully before you post any information via the Services. 6. Your Options Regarding the Information that We Collect We offer you certain choices about what information we collect from you, how we use and disclose the information, and how we communicate with you. You may choose to: Stop receiving marketing or promotional e-mails, direct mail, phone and mobile marketing communications; Update and correct the information that you provide to us; or Request that we cancel your account or that we no longer use the information you have given us to provide you services. To exercise your choices, you may: Opt-out of receiving marketing or promotional communications by following the directions in a marketing e-mail, direct mail, or mobile communication that you receive from us; Log in to your Publix Account at to update your information or change your communications preferences; Return to the area of our Services at which you signed up for a particular program and modify your information; Call 800-242-1227 with your request and current contact information; or Contact us online with your request and current contact information. If you opt-out using the directions and/or mechanisms provided in a specific type of communication, please be aware that you may only stop receiving that type of communication unless you exercise another option as described herein. Please note that if you request that your information no longer be used to provide you services, Publix or its vendors may retain your information for Publix's business purposes, such as for backups, analysis, records of deletions, legal requirements, and data retention requirements or policies. Pharmacy customers are entitled under HIPAA to exercise certain rights regarding PHI, such as access to records. For more information, please see the Pharmacy Notice Of Privacy Practices referenced above. Virginia residents have additional privacy rights under applicable law. See the Virginia Privacy Notice for additional information. 7. Cookies, Tracking Technologies and Interest-Based Advertising By using the Services, you consent to our use of cookies, web beacons, and similar technologies, including those operated by third parties on our behalf. We use these technologies to receive and store certain types of information when you interact with the Services. The information obtained by using these technologies (including, but not limited to, the pages you visit on our website, which web address you came from, information about the devices and software you use, and the content you view) helps us in various ways, including, but not limited to, to customize your experience, improve the Services, learn what areas of the Services are of most interest to our customers, and troubleshoot any issues. In order to provide the best experience possible, we also use this information for reporting and analysis purposes. Most web browsers are set to accept cookies by default, though you can manually disable the cookies feature of your browser. You should review the online documentation relating to your browser or consult with the provider of your browser for instructions on how to disable cookies. Please note that if you disable cookies on your Internet-enabled device, some parts of the Services may not function properly or offer you the same degree of personalization. In addition, you may be able to exercise choices regarding companies such as online network advertisers via certain third-party tools. Read more information about such choices at: . We may also use technologies that let us know whether you viewed certain content, whether you viewed our emails or other electronic communications, and whether the advertising we place on the Services and other websites or mobile applications is effective. One of the services we use to learn about how you and others use the Services is Google Analytics. You can learn more about how Google uses the information it collects at: . Our digital properties are not designed to respond to "do not track" requests from browsers. Visit the following website, www.allaboutdnt.org, for more information on this developing area. You can opt out of interest-based advertising from third-party providers who follow the Digital Advertising Alliance's Self-Regulatory Principles for Online Behavioral Advertising at www.aboutads.info/choices. 8. Third-party Content and Activities The Services may contain links to content or integrated content (like embedded videos) provided by third parties. The information practices of those third parties may differ from ours. You should consult the privacy policies of any third party that provides content that you access through the Services as we do not control and are not responsible for the information that they collect, use, or share. Because your web browser may obtain integrated content directly from third-party servers, the third parties may be able to collect information as if you had visited their sites directly. Some of the third parties that provide content on our Services or support our business by facilitating online advertising, analytics, or research may use technologies to track your online activities over time and across various websites and other online services. 9. Security Publix uses administrative, technical, and physical safeguards designed to safeguard the information we collect via the Services. However, no information system can be 100% secure. We cannot guarantee the absolute security of your information. Moreover, we are not responsible for the security of information you transmit or receive via the Services over networks that we do not control, including the Internet and wireless networks. You are responsible for protecting the confidentiality of your username and password(s) and you hereby release Publix from any liability related to your disclosure of your username and/or password(s) to any third party. Publix reserves the right to suspend or terminate your account or to require you to change your username and/or password if Publix believes they are no longer secure. We will make any legally required disclosures in the event of any compromise of personal information. To the extent the law allows us to provide such notification via e-mail or conspicuous posting on the Services, you agree to accept notice in that form. 10. Protecting Children's Privacy We do not knowingly solicit or collect personal information from children under the age of 13 without parental consent. If you believe that a child under age 13 may have provided us with personal information without parental consent, please contact us. 11. International Data Transfers and Use Our Services are controlled and operated by us from the United States and are not intended to subject us to the laws or jurisdiction of any state, country or territory other than those of the United States. Any information you provide through use of the Services may be stored and processed, transferred between and accessed from the United States and other countries which may not guarantee the same level of protection of personal information as the one in which you reside. However, we will handle your personal information in accordance with this Policy regardless of where your personal information is stored or accessed. 12. Contacting Us Please contact us if you have any questions about our Privacy Policy, our privacy practices, or if you would like us to update information or preferences you provided to us. Home All Products All Products DenTek® Temparin® Max Advanced Repair Kit Temporarily hold lost fillings and inlays in place until you have a chance to visit the dentist with DenTek® Temparin® Max Advanced Repair Kit. Available in two sizes: 2.64g (13+ repairs) and 1.1g (5+ repairs) Contains the same maximum-strength ingredients dentists use. Temporarily repairs lost fillings, caps, crowns, or inlays Should only be removed by a dental professional, who should be seen within 48 hours of use Directions to replace a lost filling: Wash the cavity area with warm water and do not dry. Cut with spade tool or pinch a small amount of DenTek® Advanced Repair from the vial, and roll into a ball between fingers. Slightly overfill the cavity with DenTek® Advanced Repair and tamp into place with applicator tool tip. Close your teeth together, as if you were biting down on a piece of paper, and gently chew or grind on the replaced filling a couple times to create a comfortable bite. Use a moist cotton swab to remove excess DenTek® Advanced Repair from around the area. Moisture will accelerate the set. DenTek® Advanced Repair will set softly in your mouth within a few minutes, but avoid chewing on the repair for 2-3 hours to allow the material to fully set.Directions to secure loose crowns, caps, or inlays: Try the crown, cap, or inlay first to see if you can fit it accurately back onto the tooth surface. If you are not able to do so, or if remaining tooth is cracked, then do not proceed further. If the crown has a post attached to it or if there is tooth surface remaining inside the crown, DenTek® Advanced Repair may not adhere. Carefully rinse the tooth area and the inside of the crown, cap, or inlay with warm water, and do not dry. Cut with spade tool or pinch a small amount of DenTek® Advanced Repair from the vial, roll into a ball between fingers, and place evenly inside the crown, cap, or inlay. Place the crown, cap, or inlay on the tooth and bite down, applying pressure to push into position. DenTek® Advanced Repair will begin to harden within a few minutes, but avoid chewing on the repair for at least 2-3 hours to allow the material to set. Use a toothbrush to gently clear any excess material from around the crown, cap, or inlay. Zinc Oxide Powder, Calcium Sulfate, Petroleum Jelly, Potassium alum, Paraffin, Aluminum Sulphate, Aluminum Phosphate, Menthol Crystals, Eugenol, Yellow Iron Oxide To avoid infection and/or serious injury, this product should not be used if there is throbbing pain, ulceration or swelling in the affected area. If any of these symptoms exist, consult your dentist immediately. This product should not be placed on open wounds. Do not use if allergic to Eugenol. READ ALL OF THE DIRECTIONS CAREFULLY BEFORE USE! Wash and dry your hands. Do not get unused Temparin Max wet. Use only what is needed to repair the dental emergency, and save the rest for a future application.KEEP OUT OF REACH OF CHILDREN.